



Manufactured by:
D & K Engineering
15890 Bernardo Center Drive
San Diego, CA 92127

Manufactured for:
Bioness Inc
25103 Rye Canyon Loop
Valencia, CA 91355 USA
Telephone: 800.211.9136
Email: info@bioness.com
Website: www.bioness.com

Vector Gait and Safety System™, Bioness, the Bioness Logo®, and LiveOn® are registered trademarks of Bioness Inc. in the United States or other countries | www.bioness.com

©2013 Bioness Inc.
All Rights Reserved

No part of this publication may be reproduced, transmitted, transcribed, stored in a retrieval system, or translated into any language or any computer language, in any form or by any third party, without the prior written permission of Bioness Inc.

Disclaimer

Bioness Inc and its affiliates shall not be liable for any injury or damage suffered by any person, either directly or indirectly, as a result of the unauthorized use or repair of Bioness Inc products. Bioness Inc and its affiliates do not accept any responsibility for any damage caused to its products, either directly or indirectly, as a result of use and/or repair by unauthorized personnel.




vector

BASE






Clinician's Guide



Table of Contents

List of Symbols	VII
Chapter 1: Introduction	1
Device Description	3
Chapter 2: Warnings and Cautions	3
Intended Use	3
Contraindications	3
Precautions	3
Warnings	3
Chapter 3: Patient Set-Up	5
Selecting a Harness Size	5
Donning the Patient Harness	5
Attaching the Patient Harness to the Spreader Bar	6
Chapter 4: Using the Vector Base System	9
Chapter 5: Troubleshooting and Maintenance	11
Chapter 6: Technical Specifications	13

List of Symbols

	Caution or Warning
	Consult Instructions for Use
	Re-Order Number
	Manufacturer
	Date of Manufacture

Introduction

Vector Base is a body weight support system designed to accelerate physical rehabilitation of patients with severe gait and/or balance impairment. The system unloads a fraction of the patient's weight, enabling the patient to practice walking with less than a full body weight. The patient may walk along a predetermined path or on a treadmill. As the patient walks, the Vector system remains above the patient which allows for the unload force to be applied in the vertical direction.

Device Description

The Vector Base System consists of a mechanical track, trolley with linear spring, rope, spreader bar, and patient harness. See Figure 1-1.

Mechanical Track

The mechanical track is mounted overhead, typically attached to a support beam, floor or truss, and the trolley runs on the track. The mechanical track defines the horizontal aspect of the walking path of the patient. The track can have one or more trolleys on it, it can include straight or curved portions, and it can be an open or a closed loop track.

Trolley

The trolley is attached and moves along the overhead mechanical track. The trolley has an integrated spring that allows the patient when using the system to experience a more natural gait.

Rope and Spreader Bar

The rope attaches the trolley to the spreader bar. The rope material is a 10mm diameter HMPE (High Modulus Polyethylene) and after installation is at a fixed length.

The spreader bar has two hoists with rope attached to carabiners, which attaches to the patient harness at two points. On each side of the spreader bar there are three holes built to accommodate patients of different widths, and the hoist is to be connected to one of the three points depending on the patient width.

Patient Harness

The patient harness is placed around the patient's torso, abdomen, and legs. The harness is then attached to the spreader bar.





Figure 1-1: Vector Base System.

Warnings and Cautions

Intended Use

Vector Gait and Safety System is a body weight support system designed to accelerate physical rehabilitation of patients with severe gait and/or balance impairment. The Vector Base system unloads an amount of the patient's weight thus enabling the patient to practice walking with less than a full body weight.

Contraindications

Prior to use all patients need to be deemed medically stable and appropriate by a Vector trained healthcare professional. Patients with the following diagnoses or symptoms should not use the Vector Gait and Safety System:

- Patients that weigh more than 400 lbs
- Unstable fractures
- Halo neck supports
- Uncontrolled hypertension
- Uncontrolled diabetes
- Seizures
- Severe osteoporosis

Precautions

- The spreader bar is low and can cause injury for the patient or therapist; use caution especially with tall patients.
- For correct usage, ensure that the spreader bar is not too narrow or too wide.
- Use care when assisting patients into/out of the harness as they can get tangled and fall.
- Carefully measure and fit the harness to the patient. Incorrect sizing may cause discomfort and/or injury to the patient.
- Use care when attaching the harness to the rope. If not attached correctly, there is a risk of the patient falling.

Warnings

- Do not use device on a patient that weighs more than 400 lbs.
- Do not exceed recommended weight limits or otherwise misuse the equipment as the equipment can dislocate and fall causing serious injury up to and including death.

- Do not use equipment if maintenance is required; serious injury and/or death can occur.
- Do not use equipment if the rope knot appears loose.

Patient Set-up

Selecting a Harness Size

It is important to make sure that the proper size harness is selected for a patient to ensure safety and effectiveness while using the Vector Base System. Refer to Table 3-1 for harness sizing.

Harness Size	Waist Size ¹	Harness Color Code
X-Small	24"-28"	Gray
Small	28"-32"	Yellow
Medium	32"-36"	Red
Double Medium	34"-38"	Purple
Large	36"-40"	Blue
X-Large	40"-44"	Green
2X-Large	44"-48"	White
3X-Large	48"-52"	Orange
4X-Large	52"-56"	Brown
5X-Large	56"-60"	Pink

¹ Waist Size is determined by measuring 2 inches below the umbilicus (belly button) tightly around the torso.

Table 3-1: Harness Sizing Chart.

Donning the Patient Harness

⚠ CAUTION: Inspect the harness prior to use. Do not use if harness straps are torn, if there is separation in the harness stitching, or if the Velcro is unable to adhere. Contact Bioness for a replacement harness.

1. Detach all Velcro attachments to open the vest.
2. Make sure the patient is in a secured position.
3. Have the patient put their arms through the shoulder straps to put on the vest. The open side of the vest will face the front.
4. Secure the lower abdominal support first by attaching the Velcro. Make sure the support is as tight as possible.
5. Next secure the upper abdominal support by attaching the Velcro. Note: This is the middle support with the netting material on the back.
6. Securely attach the top torso support.
7. Tighten side straps on lower abdominal support.
8. If necessary adjust side vertical straps.
9. Wrap leg supports around the upper thigh and secure the Velcro.

10. If necessary adjust side vertical straps attached to the leg supports.
11. Adjust shoulder straps to the desired length.




Figure 3-1: Patient Harness Features.

Attaching the Patient Harness to the Spreader Bar

After the patient harness has been securely placed and adjusted on the patient, the harness will need to be attached to the spreader bar.

1. The spreader bar provides three attachment points to attach the harness. The selection of the attachment point depends on the width of the patient. Connect the rope hoist to the outermost attachment points for wide patients, and to the innermost points for narrow patients. See Figure 3-2.
2. Release the rope hoist lever to lower the side of the rope with the attached carabiner.
3. Connect the carabiner to the clip on the patient harness shoulder strap.
4. Pull the side of the rope that is not connected to the carabiner until the rope connected to the harness shoulder strap is tight.
5. Repeat steps 1-3 for the other patient harness shoulder strap.
6. Pull on the ropes that are not connected to the harness shoulder straps to adjust the level of tension in the harness shoulder straps.
7. After the shoulder straps have been properly adjusted, tuck the loose ropes into the back of the harness.

Note: While using the Vector Base System the harness may slide up on the patient causing a decrease in rope tension. If this happens un-tuck the loose ropes from the back of the harness, adjust the tension level, and tuck the loose ropes into the back of the harness.




Figure 3-2: Spreader Bar attachment points.

Using the Vector Base System

Safety Features

The Vector Base System includes the following safety features:

- Mechanical tracks that are in an open configuration contain stop bumpers to prevent the trolley from coming off the end of the track.
- If a patient falls while attached to the Vector Base System, the trolley has been designed to provide cushioning so that the patient does not experience a severe jerk at the end of the rope.

Rehabilitation

In a rehabilitation environment the Vector Base System provides body weight support to patients for safe gait training, balance exercises and treadmill training. The system is designed to hold a patient up to 400 lbs and provide up to 180 lbs of body weight support via the use of a linear spring, which is built into the trolley.

The patient's path of travel is limited to the path of an overhead mechanical track which the trolley is attached to. The system is designed to provide a safe rehabilitation environment for both patients and healthcare providers through the continuum of care.

Troubleshooting and Maintenance

If you have any questions or concerns about the Vector Base System, please contact Bioness at (800) 211-9136, Option 3.

Troubleshooting

Issue	Problem
Trolley Squeaking	Inspect track and drive wheels for debris and dirt build up. Wipe off with a dry cloth.
Harness Will Not Fit Patient Properly	<ul style="list-style-type: none"> • Verify that appropriate size is being used. • If appropriate use a larger harness. • Inspect Velcro and material to ensure integrity is still intact.
Rope is Worn, Frayed or Damaged	<ul style="list-style-type: none"> • Replace rope. • Call Bioness to order replacement part.

Maintenance

Routine maintenance is critical for the safety and effectiveness of this device.

Check the following items before each use of the Vector Base System:

- Inspect the rope. If the rope is worn, frayed or damaged replace before use.
- Inspect the patient harness for tears, material separation and Velcro degradation.
- Inspect the spreader bar before each use.
- Inspect the carabiners before each use.
- Inspect the mechanical track quarterly. Clean debris and/or dirt when necessary. Tighten loose bolts, screws or hardware if appropriate.

Items to be inspected and maintained on an annual basis:

- Visually inspect the trolley wheels for wear.
- Measure that the trolley spring is fully compressed when loaded with 180 lbs.
- Replace the rope connecting the trolley to the spreader bar annually. Call Bioness at (800) 211-9136 to order a replacement part.
- Replace the rope hoists with attached carabiners annually. Call Bioness at (800) 211-9136 to order a replacement part.

Cleaning Instructions for the Patient Harness

1. Hand wash in cold water.
2. Hang dry (do not place harness in dryer).
3. Do not iron.

Technical Specifications

Vector Base System Technical Specifications	
Trolley Dimensions	8.75" long, 4.5 inches wide, and 10.4 inches tall
Trolley Material	Aluminum grade 6061-T6
Trolley Patient Weight Limit	400 lbs
Rope Material	10mm diameter HMPE (High Modulus Polyethylene)
Patient Harness Material	Lycra, Rayon, Cotton and Neoprene